Levenhuk 320\320L Biological Microscopes

User Manual

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levenhuk
Zoom&Joy

1 Eyepiece
 2 Monocular head
 3 Revolving nosepiece
 4 Arm
 5 Objective lenses
 6 Specimen holder
 7 Stage
 8 Focus adjustment knob
 9 Condenser
 10 Coaxial stage control
 11 Illumination
 12 Base

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Never look directly at the Sun through optical devices or with the unaided eye! It can result in permanent eye damage or even BLINDNESS.

General use

Levenhuk 320\320L biological microscopes are designed for observations of specimens in transmitted light. They may be used in many scientific fields, including biology, ecology, medicine and veterinary science, as well as for technological research and school presentations.

Parts of the microscope

Base. It supports the weight of the microscope and houses the illumination source, electronics and control mechanisms.

Arm. This piece holds the base, the stage and the head of the microscope together. This microscope has coaxial coarse and fine focusing systems with an adjustment lever and a locking mechanism, which provides for smooth vertical movements of the stage.

Rack-and-pinion mechanism. The stage with the condenser is moving along this column, which is mounted on the arm of the microscope. For additional precision, a condenser may be aligned separately.

Eyepieces. Wide FOV eyepieces (WF10x and WF16x) are used in this microscope.

Revolving nosepiece. This quadruple nosepiece allows you to change objective lenses smoothly and easily.

Objective lenses. High-quality achromatic objective lenses with 4x, 10x, 40x (S) and 100x (S, oil immersion) magnifications provide for sharp and bright images.

Stage. Sturdy and reliable stage with a specimen holder that can be used to move your slides while observing them.

Condenser. Abbe condenser with 1.25 N.A. iris diaphragm.

Using the microscope

- Insert an eyepiece into the monocular head and attach objective lenses to the revolving nosepiece in ascending order. Place a specimen in the middle of the stage and hold it in place with stage clips.
- Connect your microscope to a power outlet and adjust the brightness of illumination.
- When you finish your observations, always dim the illumination before unplugging the instrument.
- Start at the lowest magnification and center the observed specimen. Afterwards, switch to a higher magnification by rotating the revolving nosepiece. If the image is out of focus, you can use the fine focus system. When using a 100x objective lens, you have to fill the space between the lens and the specimen with immersion oil. There should be no air bubbles in the oil, as this will reduce the resulting quality of the image. For the best quality of your views, illumination has to be carefully adjusted. When switching to a different objective lens, you will most likely have to readjust the iris diaphragm and the illumination brightness.
- If you need to change the light bulb, unplug your microscope, wait for the bulb to cool and change it. Note that the lamp has to be leveled and inserted all the way to the end.

Using the oil immersion objective lens

Before observing at high magnifications with a 100x objective lens, you have to place a drop of immersion oil on the specimen and the objective lens itself. You can do it with a glass tube or a dropper. A vial of oil is supplied with the microscope.

Do not use any substitutes for the immersion oil, as this will significantly reduce the quality of the resulting image. After placing the oil on the lens and the specimen, raise the microscope stage with a coarse focus knob until two drops merge. To bring the image into focus, you can use the fine focus system.

When you are finished observing, clean the immersion oil with a cloth or a filter paper. If the objective lens has become stained you can clean it with a cotton pad dipped in alcohol. **Be very gentle and do not apply any excessive pressure, as you may damage the objective lens!**

Specifications

Magnification 40x-1600x
Objective lenses 4x, 10x, 40x and 100x
Revolving nosepiece quadruple
Eyepieces WF10x, WF16x
Condenser Abbe, with 1.25 N.A. iris diaphragm
Illumination 12W halogen lamp
Focus system coaxial, coarse and fine

The kit

- 4x, 10x, 40x, 100x (Oil) objective lenses Wf10x
- Wf16x eyepieces
- Abbe condenser
- Immersion oil
- D320L: Levenhuk C310 NG 3 mpx digital camera, MiniSee ScopePhoto installation CD-ROM
- User manual and warranty

Digital camera

Levenhuk D320L comes with a 3 mpx C310T NG digital camera.

The camera allows you to observe specimens in fine detail and true colors on your PC monitor and save images on the hard drive.

The special software that comes in the kit allows you to view and edit the resulting images. Supported file formats include: *.bmp, *.jpg, *.jpeg, *.png, *.tif, *.tiff, *.gif, *.psd, *.ico, *.emf, *.wmf, etc.

Connect the camera to your PC via a USB-cable that also acts as a power cable, so you do not have to plug it into a power source.

Compatible OS: Windows 2000/XP/2003/Vista

Camera model C310 NG
Max. resolution (for still images) 2048x1536
Megapixels 3.0
Sensor 1/2° CMOS
Dimensions (cylinder unit maximum diameter), mm 60
Pixel size 3.2µm x 3.2µm
Sensitivity, V/lux.sec@550 nm 1.5
Mounting location eyepiece tube (replaces the eyepiece)
Image format BMP, TIFF, JPG, PICT, SFTL, etc.
Field of view diameter 18 mm
Spectral range 400 nm - 650 nm
Exposure ERS
White balance auto/manual
Exposure control auto/manual
Software image control image size, brightness, exposure control
Interface USB 2.0, 480 Mb/s
System requirements Windows 2000/XP(SP2)/2003/Vista; USB port

Software USB 2.0 driver, software for image viewing and editing
Shell Metal cylinder, painted black
Power source USB 2.0 cable
Temperature range -22... 158 °F (-30... 70 °C)

Caution

Never, under any circumstance, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS. Take necessary precautions when using the device with children or people, who have not read and fully understood these instructions.
Do not try to take the device apart on your own. For repairs of any kind, please contact your local specialized service center.

Batteries safety instructions

- Always purchase the correct size and grade of battery most suitable for the intended use.
- Always replace the whole set of batteries at one time; taking care not to mix old and new ones, or batteries of different types.
- Clean the battery contacts and also those of the device prior to battery installation.
- Make sure the batteries are installed correctly with regard to polarity (+ and -).
- Remove batteries from equipment that is not to be used for an extended period of time.
- Remove used batteries promptly.
- Never attempt to recharge primary batteries as this may cause leakage, fire, or explosion.
- Never short-circuit batteries as this may lead to high temperatures, leakage, or explosion.
- Never heat batteries in order to revive them.
- Remember to switch off devices after use.
- Keep batteries out of the reach of children, to avoid risk of ingestion, suffocation, or poisoning.

Care and maintenance

- After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection.
- Do not touch the optical surfaces with your fingers. You can use a soft napkin or a cleaning wipe, dipped in ether or absolute alcohol, to clean lenses. We recommend using special optics cleaning fluids from Levenhuk.
- Do not use any corrosive or acetone-based fluids to clean the optics.
- Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush.
- Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws. Protect the device from sudden impacts and excessive mechanical force.
- Store the device in a dry, cool place away from dust, hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures. Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity.
- Be careful during your observations, always replace the dust cover after you are finished with observations to protect it from dust and stains.
- If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope.
- Seek medical advice immediately if a small part or a battery is swallowed.

Levenhuk Limited Warranty

All Levenhuk telescopes, cameras for telescopes, microscopes and binoculars are warranted to be free of defects in materials and workmanship for **three years** from date of retail purchase. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from date of retail purchase. Levenhuk will repair or replace such product or part thereof which, upon inspection by Levenhuk, is found to be defective in materials or workmanship. As a condition to the obligation of Levenhuk to repair or replace such product, the product must be returned to Levenhuk together with proof of purchase satisfactory to Levenhuk.

This warranty does not cover consumable parts, such as batteries.

A Return Authorization (RA) Number must be obtained in advance of return. Contact the local Levenhuk branch to receive the RA number to be displayed on the outside of your shipping container. All returns must be accompanied by a written statement setting forth the name, address and telephone number of the owner, including a description of any claimed defects. Parts or products for which replacement is made will become the property of Levenhuk.

The customer will be responsible for all costs of transportation and insurance to and from Levenhuk or its authorized dealers and will be required to prepay such costs.

Levenhuk will use reasonable efforts to repair or replace any product covered by this warranty within thirty days of receipt. If a repair or replacement will require more than thirty days, Levenhuk will notify the customer accordingly. Levenhuk reserves the right to replace any product that has been discontinued from its product line with a new product of comparable value and function.

This warranty does not apply to any defects or damages resulting from alteration, modification, neglect, misuse, usage of improper power sources, damage in transportation, abuse, or any cause other than normal use, or to malfunction or deterioration due to normal wear.

Levenhuk disclaims all warranties, express or implied, whether of merchantability or fitness for a particular use, except as expressly set forth herein. The sole obligation of Levenhuk under this limited warranty will be to repair or replace the covered product, in accordance with the terms set forth herein. Levenhuk disclaims liability for any loss of profits, loss of information, or for any general, special, direct, indirect or consequential damages which may result from breach of any warranty, or arising out of the use or inability to use any Levenhuk product. Any warranties which are implied and which cannot be disclaimed will be limited in duration to a term of three years for telescopes, cameras for telescopes, microscopes and binoculars or six months for accessories from the date of retail purchase.

Some states/provinces do not allow the exclusion or limitation of incidental or consequential damages, so the above limitations and exclusions may not apply to you. This warranty gives you specific legal rights, and you may have other rights which vary from state to state or province to province.

Levenhuk reserves the right to modify or discontinue any product without prior notice.

NOTE: This warranty is valid to USA and Canadian customers who have purchased this product from an authorized Levenhuk dealer in the USA or Canada. Warranty outside the USA or Canada is valid only to customers who purchased from an authorized Levenhuk dealer in the specific country or international distributor. Please contact them for any warranty service.

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch:

Levenhuk Worldwide: USA: www.levenhuk.com, Czech Republic: www.levenhuk.cz, Russia: www.levenhuk.ru, Ukraine: www.levenhuk.com.ua, EU: www.levenhuk.eu

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Levenhuk® is a registered trademark of Levenhuk, Ltd.
Copyright © 2006-2013 Levenhuk Ltd. All rights reserved
Developed and manufactured for Levenhuk, Ltd., IL 60139 USA.