Levenhuk D70L Digital biological microscope

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levenhuk
Zoom&Joy

- 1 Monitor
- 2 LCD screen
- 3 Tubus
- 4 Revolving nosepiece
- 5 Objective lenses
- 6 Arm
- 7 Stage
- 8 Diaphragm disc
- 9 Illumination
- 10 Focus adjustment knob
- 11 Illumination adjustment wheel
- 12 Base
- 13 AC jack

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Levenhuk D70L microscope

General Information

This digital microscope is designed for bright field observations of transparent and solid objects in either transmitted or reflected light. It is well-suited for school studies in biology and may be used for microphotography.

The optics yields sharp and bright views throughout the field of view. The standard kit of this microscope even allows you to edit the microphotography images that you take during research. This microscope is manufactured in accordance with all existing regulations and technical documentation; and was certified to be ready for personal and professional use.

Specifications

Magnification	40x-400x
Objective lenses	4x, 10x, 40x
Linear field	18 mm
Tubus length	160 mm (6.3 in)
Stage size	95 mm x 95 mm (3.7 x 3.7 in)
Stage movement range	0-15 mm (0-0.6 in)
Illumination source power supply	Three AA batteries or power cable
Illumination source (lower/upper)	LED lamp
LCD monitor	3.6"
Matrix	2 Mpx CMOS (1600x1200, 800x600)
Image format	JPG, AVI
Interface	USB, AV

Package

- 4x, 10x, 40x objective lenses
- Diaphragm disc
- LCD monitor with SD card slot
- Microtome
- Tweezers
- Brine shrimp hatchery
- Flasks with yeast, brine shrimp eggs, pitch, and sea salt
- Prepared samples (linden, maize, onion, pine, pumpkin, bean, tree bark, cotton, mosquito)

Parts of the microscope

Base: supports the weight of the microscope and houses the illumination source, electronics and control mechanisms.

Arm: holds the base, the stage and the head of the microscope together.

Rack-and-pinion mechanism: is mounted on the microscope arm, provides for vertical movements of the stage. Coaxial control knobs are placed on both sides of the arm. Total adjustment range is 15 mm.

Stage: is equipped with a diaphragm disc and a specimen holder that can be used to move slides during your observations.

Locking screw: to avoid any damage to the objective lens, the microscope arm has a locking screw that acts as an upper limit for stage movements.

Revolving nosepiece: houses three objective lenses of the microscope, and allows you to change these lenses smoothly and easily.

LCD monitor: instead of an eyepiece this microscope uses a 3.6" LCD monitor.

Choosing the correct objective lens

Always begin your observations at the lowest magnification possible. This will allow you to center the specimen in your field of view. This alignment of the specimen has to be done carefully, or you might not get the desired view upon changing the magnification.

Focusing the view

- Turn on the LCD monitor.
- Place a specimen on the stage and use stage clips to hold it in place.
- Turn the revolving nosepiece until 4x objective lens is locked in place. Move the slide manually, to center the desired part of the observed object in the field of view.
- Raise the stage by turning the focus knobs until the cover glass of the slide touches the objective lens. Be careful and try not to damage the slide or the lens.
- Look at the monitor and, if necessary, adjust focus with a fine focus knob until you see a clear image on the screen.

If this procedure is performed properly, you will not have to worry about damaging the optics or the specimen, when switching to a different objective lens. However, additional focusing may still be required.

Environment conditions

This model performs best under the following conditions: room temperatures from 10 °C to 35 °C (50 °F to 95 °F), relative humidity not exceeding 80% at 25 °C (77 °F).

The room should also be void of any dust, acid or base vapors and other active substances

Caution

Never, under any circumstance, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS.

Take necessary precautions when using the device with children or people, who have not read and fully understood these instructions.

Do not try to take the device apart on your own. For repairs of any kind, please contact your local specialized service center.

Care and maintenance

- After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection.
- Do not touch the optical surfaces with your fingers. You can use a soft napkin or a cleaning wipe, dipped in ether or absolute alcohol, to clean lenses. We recommend using special optics cleaning fluids from Levenhuk. Do not use any corrosive or acetone-based fluids to clean the optics.
- Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush.
- Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws. Protect the device from sudden impacts and excessive mechanical force.

- Store the device in a dry, cool place away from dust, hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures. Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity.
- Be careful during your observations, always replace the dust cover after you are finished with observations to protect it from dust and stains.
- If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope.

Troubleshooting

If you are experiencing some problems with your instrument during its service, refer to this section for possible solutions. However, if this doesn't help, you should contact your local service center.

Problem	Possible cause	Solution
Uneven illumination throughout the field of view	The objective lens is not locked in place, i.e. optical system is not properly aligned.	Rotate the revolving nosepiece until the lens locks in place.
	The diaphragm disc is not locked in place	Rotate the disc until it locks in place
	One of the lenses is dirty	Clean the lens from dust or dirt
Clean the stage from dust	One of the lenses is dirty	Clean the lens from dust or dirt
	The stage is dirty	Clean the stage from dust
Reduced quality of views (low resolution, poor contrast)	One of the lenses is dirty	Clean the lens from dust or dirt
The field of view is clipped	The objective lens is not locked in place	Rotate the revolving nosepiece until the lens locks in place
	The diaphragm disc is not locked in place	Rotate the disc to lock the diaphragm in place
	The slide is moving	Use stage clips to hold it in place
When changing objective lenses, they scrape the slide on the stage	The slide is upside down	Place it correctly
	The cover glass is too thick	Use a cover glass of standard thickness
The illumination system does not work	One of the batteries is needs replacing	One of the batteries is needs replacing Change
The monitor does not work	The power cable is not plugged in	Make sure that the power cable is plugged into the socket and the AC jack

International Warranty

- The Vendor guarantees that the quality of the Levenhuk product purchased complies with the technical documentation requirements on conditions that the consumer observes rules of transportation, storage conditions and operating instructions.
- Levenhuk Ltd. warrants the product against defects in materials.
- During the warranty period the Buyer can return the defect product to the Vendor or to Levenhuk service center. Levenhuk Ltd. will repair or replace the defect product at its own discretion.
- No claims are accepted in case the properly filled warranty slip is absent or contains corrections, or if the defect product has not been provided by the Buyer. The manufacturer or the seller is not liable for any damages caused by misuse of the product. No guarantee is provided in cases the product was used for purposes other than that intended or has mechanical damages, scratches, cracks, optics damages, as well as in cases the product is broken-down due to impact damages, squeezing or stretching or the product has been disassembled or repaired by unauthorized personnel.
- This warranty does not cover components and/or accessories with limited use period, batteries, etc.

The period of warranty is 3 (three) years beginning on the date of purchase. Please keep the warranty slip along with your receipt.

For any warranty services please contact the local Levenhuk branch.

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Levenhuk® is a registered trademark of Levenhuk, Ltd.
Copyright © 2006-2013 Levenhuk Ltd. All rights reserved
Developed and manufactured for Levenhuk,Ltd.